




Hydrophilic Guidewire

国械注准20213030017



Application

The hydrophilic coated guidewire with special coating formula is widely used in both diagnosis and treatment of interventional therapies.

It has good support yet right control to ensure device delivery through tortuous anatomy.


Features

- Advanced hydrophilic coating to provide lasting lubrication.
- Soft atraumatic tip provides flexibility and minimizes trauma.
- The NT core provides 1:1 torque control, perfect combination of support and maneuverability in tortuous anatomy.
- Enhanced visualization for accurate positioning.
- Complete range of specification to accommodate physician preference.

Ordering Information

NO.	Model NO.	Diameter (inch)	Length (cm)	Tip configuration
1	GWO - 14 - 150S	0.014	150	Straight
2	GWO - 14 - 150A		150	Angled
3	GWO - 18 - 120A	0.018	120	Angled
4	GWO - 18 - 150S		150	Straight
5	GWO - 18 - 150A		150	Angled
6	GWO - 18 - 180S		180	Straight
7	GWO - 18 - 180A		180	Angled
8	GWO - 18 - 220S		220	Straight
9	GWO - 18 - 220A		220	Angled
10	GWO - 18 - 260S		260	Straight
11	GWO - 18 - 260A		260	Angled
12	GWO - 21 - 120A	0.021	120	Angled
13	GWO - 21 - 150S		150	Straight
14	GWO - 21 - 150A		150	Angled
15	GWO - 21 - 180S		180	Straight
16	GWO - 21 - 180A		180	Angled
17	GWO - 21 - 220S		220	Straight
18	GWO - 21 - 220A		220	Angled
19	GWO - 21 - 260S		260	Straight
20	GWO - 21 - 260A		260	Angled
21	GWO - 25 - 120A	0.025	120	Angled
22	GWO - 25 - 150S		150	Straight
23	GWO - 25 - 150A		150	Angled
24	GWO - 25 - 180S		180	Straight
25	GWO - 25 - 180A		180	Angled
26	GWO - 25 - 220S		220	Straight
27	GWO - 25 - 220A		220	Angled
28	GWO - 25 - 260S		260	Straight
29	GWO - 25 - 260A		260	Angled
30	GWO - 32 - 120A	0.032	120	Angled
31	GWO - 32 - 150S		150	Straight
32	GWO - 32 - 150A		150	Angled
33	GWO - 32 - 180S		180	Straight
34	GWO - 32 - 180A		180	Angled
35	GWO - 32 - 220S		220	Straight
36	GWO - 32 - 220A		220	Angled
37	GWO - 32 - 260S		260	Straight
38	GWO - 32 - 260A		260	Angled
39	GWO - 35 - 050A	0.035	50	Angled
40	GWO - 35 - 080A		80	Angled
41	GWO - 35 - 120A		120	Angled
42	GWO - 35 - 150S		150	Straight
43	GWO - 35 - 150A		150	Angled
44	GWO - 35 - 180S		180	Straight
45	GWO - 35 - 180A		180	Angled
46	GWO - 35 - 220S		220	Straight
47	GWO - 35 - 220A		220	Angled
48	GWO - 35 - 260S		260	Straight
49	GWO - 35 - 260A		260	Angled
50	GWO - 38 - 050A	0.038	50	Angled
51	GWO - 38 - 080A		80	Angled
52	GWO - 38 - 120A		120	Angled
53	GWO - 38 - 150S		150	Straight
54	GWO - 38 - 150A		150	Angled
55	GWO - 38 - 180S		180	Straight
56	GWO - 38 - 180A		180	Angled
57	GWO - 38 - 220S		220	Straight
58	GWO - 38 - 220A		220	Angled
59	GWO - 38 - 260S		260	Straight
60	GWO - 38 - 260A		260	Angled

Models Naming Composition



E.g. Model No.GWO-35-150A represents guidewire OD of 0.035 inch, 150cm in length, with angled tip.

- Commonly used models, recommend to stock up regularly.
- Less commonly used models, recommend to stock up moderately.

Note: 0.01 inch = 0.254 mm