

# Angiographic Catheter



**国械注准20213030823**

## Application

Angiographic Catheters are designed to provide a pathway for delivering contrast media to selected sites in the vascular system.

## Features

- Soft tip material designed to reduce risk of trauma during contrast injection.
- Highly radiopaque material promotes clear flourosopic visualization.
- Durable, fluted hub designed to ease guidewire manipulation and provide strain relief for kink-resistant advancement.
- Wide variety of shapes provides the versatility required to meet individual clinical requirements.

Ordering Information

| Terminal configuration     | Shapes                                                                              | Models         | OD (Fr) | ID (mm) | length (cm) |
|----------------------------|-------------------------------------------------------------------------------------|----------------|---------|---------|-------------|
| Hepatic artery tube (RH)   |    | RH630-4.2-80   | 4.2     | 1.0     | 80          |
|                            |                                                                                     | RH630-5-80     | 5       | 1.03    | 80          |
| Simmons tube (SM)          |    | SM510-4.2-100  | 4.2     | 1.0     | 100         |
|                            |                                                                                     | SM510-5-100    | 5       | 1.03    | 100         |
|                            |    | SM511-4.2-100  | 4.2     | 1.0     | 100         |
|                            |                                                                                     | SM511-5-100    | 5       | 1.03    | 100         |
|                            |    | SM512-4.2-100  | 4.2     | 1.0     | 100         |
|                            |                                                                                     | SM512-5-100    | 5       | 1.03    | 100         |
| Hunter head tube (HH)      |    | HH371-4.2-100  | 4.2     | 1.0     | 100         |
|                            |                                                                                     | HH371-5-100    | 5       | 1.03    | 100         |
| Renal artery tube (HS)     |  | HS640-4.2-80   | 4.2     | 1.0     | 80          |
|                            |                                                                                     | HS640-5-80     | 5       | 1.03    | 80          |
|                            |                                                                                     | HS641-4.2-80   | 4.2     | 1.0     | 80          |
|                            |                                                                                     | HS641-5-80     | 5       | 1.03    | 80          |
| Cobra tube (CB)            |  | CB230-4.2-80   | 4.2     | 1.0     | 80          |
|                            |                                                                                     | CB230-5-80     | 5       | 1.03    | 80          |
|                            |                                                                                     | CB231-4.2-80   | 4.2     | 1.0     | 80          |
|                            |                                                                                     | CB231-5-80     | 5       | 1.03    | 80          |
|                            |                                                                                     | CB232-4.2-80   | 4.2     | 1.0     | 80          |
|                            |                                                                                     | CB232-5-80     | 5       | 1.03    | 80          |
| Single bend tube (RC)      |  | RC620-4.2-80   | 4.2     | 1.0     | 80          |
|                            |                                                                                     | RC620-5-80     | 5       | 1.03    | 80          |
|                            |                                                                                     | RC621-4.2-80   | 4.2     | 1.0     | 80          |
|                            |                                                                                     | RC621-5-80     | 5       | 1.03    | 80          |
|                            |                                                                                     | RC622-4.2-80   | 4.2     | 1.0     | 80          |
|                            |                                                                                     | RC622-5-80     | 5       | 1.03    | 80          |
| Multifunctional tube (MPA) |  | MPA35-5-110    | 5       | 1.03    | 100         |
|                            |                                                                                     | MPA280-4.2-100 | 4.2     | 1.0     | 100         |
|                            |                                                                                     | MPA280-5-100   | 5       | 1.03    | 100         |
|                            |                                                                                     | MPA281-4.2-100 | 4.2     | 1.0     | 100         |
| vertebral artery tube (AN) |  | AN-4.2-100     | 4.2     | 1.0     | 100         |
|                            |                                                                                     | AN-5-100       | 5       | 1.03    | 100         |
| Splenic artery tube (RS)   |  | RS631-4.2-80   | 4.2     | 1.0     | 80          |
|                            |                                                                                     | RS631-5-80     | 5       | 1.03    | 80          |
| Shepherd tube (SH)         |  | SH210-4.2-80   | 4.2     | 1.0     | 80          |
|                            |                                                                                     | SH210-5-80     | 5       | 1.03    | 80          |
|                            |                                                                                     | SH211-4.2-80   | 4.2     | 1.0     | 80          |
|                            |                                                                                     | SH211-5-80     | 5       | 1.03    | 80          |
|                            |                                                                                     | SH212-4.2-80   | 4.2     | 1.0     | 80          |
|                            |                                                                                     | SH212-5-80     | 5       | 1.03    | 80          |

Models Naming Composition

XXXX-XX-XXX



E.g. Model No. SH210-4.2-80 represents catheter type of Sheperd210, ,OD of 4F, 80 cm in length.

Note: 1F = 0.33 mm